


## Biacore T200 Single-Cycle Kinetics Protocol

**STOP**

This protocol is for qualified users operating Biacore T200 in IBC 402 only. Dr. Jao accepts no responsibility for actions taken as a result of using this protocol. Reading the manufacturer's handbooks is highly recommended.

1. In this protocol, we will use the Biacore single-cycle kinetics (SCK) method to obtain kinetics data, using samples from Biacore T200 Getting Started kit. Choose **Run:Method**. In the folder **Methods And Templates**, select **Biacore Methods**. Double click **Single-cycle kinetics** to open the default setting for SCK.
2. Check the **Overview**. The overview shows a summary of the method.
3. Click the **General Setting** and set the parameters as followed:  
Data collection rate: 10Hz,  
Detection mode: Dual (choose "Multiple" if all 4 channels are needed),  
Temperature: 25,  
Concentration unit: nM  
Buffer A: HBS-EP+
4. Click the **Assay Steps** and select the steps to alter the setting.  
Select the **Startup** and change the **Number of replicates** to 3. No changes for the **Sample** step.
5. Click the **Cycle Types** and select the **Sample** cycle. The cycle is executed sequentially according to the commands. Under **Commands** panel, select **Sample 1** and set the parameters as followed:  
Type: Single cycle kinetics  
Concentration per cycle: 5  
Contact time: 120 s  
Dissociation time: 600 s  
Flow rate: 30  $\mu$ L /min  
Flow path: Both  
  
Under **Commands** panel, select the **Regeneration 1** and set the parameters as followed:  
Regeneration solution: 10mM Glycine pH2.5  
Contact time: 30s  
Flow Rate: 30  $\mu$ L /min  
Flow path: Both  
Stabilization period: 30s



6. Still under the **Cycle Types**, select the **Startup** cycle. Under **Commands** panel, select **Sample 1** and set the parameters as followed:

Type: Low sample consumption

Contact time: 120 s

Dissociation time: 120 s

Flow rate: 30 µL /min

Flow path: Both

Under **Commands** panel, select the **Regeneration 1** and set the parameters as followed:

Regeneration solution: 10mM Glycine pH2.5

Contact time: 30s

Flow Rate: 30 µL/min

Flow path: Both

Stabilization period: 30s

7. Click on **Variable Settings**. Check  Define all values at run time for both Startup and Sample steps.
8. Click on **Verifications** to verify the method.
9. Click on **Setup Run**. Choose **Flow path**. Click on **Next**. Use HBS-EP+ as **Startup** solution. Select **Sample** and type in *beta2-microglobulin* as the sample solution and the concentrations are 0 nM, 0 nM, 0 nM, 0 nM and 0 nM in the first and second cycles and 2 nM, 4 nM, 8 nM, 16 nM and 32 nM in the third. The **MW** of beta2-microglobulin is *11800* Da. Click on **Next** when done.
10. The **Cycle run list** shows a summary of the experiment cycles. Click on **Next**.

11. In **System Preparations**, uncheck **Prime before run** as well as **Normalize detector**. Click on **Next**.
12. Prepare your samples: Dilute stock beta2-microglobulin in running buffer to 32 nM (2.26  $\mu$ L beta2-microglobulin + 598  $\mu$ L HBS-EP). Prepare the concentration series from the 32 nM sample: mix 300  $\mu$ L of the 32 nM solution with 300  $\mu$ L running buffer to get the 16 nM solution. Continue the dilution series to obtain the following: 32, 16, 8, 4 and 2 nM.

Position	Volume ( $\mu$ l)	Content	Type	Sample 1 Conc (nM)	Sample 1 MW (Da)
R2 F1	320	Glycine pH2.5	Regeneration		
R2 E1	154	HBS-EP+	Startup		
R2 C4	118	beta2-microglobulin	Sample	32	11800
R2 C3	118	beta2-microglobulin	Sample	16	11800
R2 C2	118	beta2-microglobulin	Sample	8	11800
R2 C1	118	beta2-microglobulin	Sample	4	11800
R2 B6	118	beta2-microglobulin	Sample	2	11800
R2 B5	118	beta2-microglobulin	Sample	0	11800
R2 B4	118	beta2-microglobulin	Sample	0	11800
R2 B3	118	beta2-microglobulin	Sample	0	11800
R2 B2	118	beta2-microglobulin	Sample	0	11800
R2 B1	118	beta2-microglobulin	Sample	0	11800

13. Prepare and position samples according to **Rack Positions**. Click on **Next**.
14. Make sure everything is correct according to the **Prepare Run Protocol** and click **Start** to begin the experiment.
15. Enter a file name for the resulting sensorgram.
16. Evaluate result using Biacore T200 evaluation software.

Results of 2009/02/25:

	Multi-Cycle	Single-Cycle
$k_a$ (1/Ms)	1.03E+6	9.10E+5
$k_d$ (1/s)	0.0030	0.0028
$K_D$ (1/M)	3.0E-9	3.1E-9

