Simplified Refeyn TwoMP protocol



This protocol is for qualified users operating TwoMP in IBC 402 only. Dr. Jao accepts no responsibility for actions taken as a result of using this protocol. Reading the manufacturer's User's Manual is highly recommended.

- I. Instrument Warm Up
 - A. Set air conditioner to 21°C and keep the door open.
 - B. Empty the water bucket inside the dehumidifier.
 - C. Open the lid on the TwoMP. Ensure that the objective and sample carrier are clean.
 - D. Turn ON the vibration-isolation table by pressing the buttons "Power" and then "Isolation".
 - E. Turn ON the switch that is located at the back of the Electronics Unit.
 - F. Turn ON the PC-1 and start "AcquireMP" **1** at least one hour before starting your experiments
 - G. Select "Regular mode" \rightarrow Select "Desired image size (Small/Regular/Large)".

5			Measurement mode	Image size	
Open Project Folder	Concentration Calculator	Replay Movie	Normal 👻	Small	•

- H. Click "Open Project Folder", set Desktop:/Temporary Data as your project folder.
- II. Preparing and Mounting Coverslips
 - A. Clean both sides of coverslip by "ddH2O \rightarrow isopropanol \rightarrow ddH2O \rightarrow isopropanol \rightarrow ddH2O" and dry it with compressed N2. Pre-cleaned coverslip can be stored in coverslip box for a short time.
 - B. Using soft-tipped tweezers, place a gasket in the middle of the aliment tool, and then place the precleaned coverslip slide on top.
 - C. Press gently on the coverslip between the wells with soft-tipped tweezers.
 - D. Open the lid on the TwoMP unit and place a small drop of immersion oil on the objective lens.
 - E. Flip and mount the prepared coverslip (with the gasket facing the top), and place two magnets diagonally at each end of the coverslip.
 - F. Select a well (i.e. upper left corner on the gasket) that you want to start your experiment on. Use lateral control at <u>speed 50x</u> to move the stage so that the red laser beam is focused on the middle of the selected well. Once this is done, you may now <u>set the well position</u> on the software. This will allow for automatic stage positioning to another well when you run your next sample.



III. Data Collection

Before you start your experiments, please ensure that your buffer, samples and calibration standards are ready at room temperature. Concentration calculator in AquireMP is helpful when you need to dilute your samples. The first sample you measured should be a calibration standard, appropriate to your experiment type (measure Calibrants in your sample's buffer).

- A. Open the lid and add 18 µl of buffer to the selected well.
- B. Lower the lid of TwoMP and click "Droplet-dilution" to find focus.



Switch between native and ratiometric views to assess image quality. C.





- D. Examine the autofocus ring. It should be a thin, clean outline with no interruptions.
- E. Inject 2 µl of your diluted sample (50 nM-200 nM) into the buffer droplet used to find focus. Mix thoroughly by pipetting up and down.
- F. Close the lid. Wait for the signal to stabilize. Press "Record" button. The software will now record

data for 1min. Do not touch or shake the table and the instrument during recording.

- G. Click "Next well" to automatically move to next well.
- H. Repeat finding focus and collecting sample data.
- IV. Data Analysis

AquireMP and DiscoverMP cannot be opened at the same time. Analysis your data on "PC-2".

A. Turn on PC-2 and ensure "Warmhole switch"



- B. Copy your files from "PC-1" and paste to your folder in "PC-2".
- C. Initiate the DiscoverMP software.
- D. Upload your data files by click "+" button or directly drag them into data file zone.
- E. Adjust the contrast limits on a particular region of the data.
 - 1. Set the upper contrast to 0.
 - 2. Shorten the lower contrast limit to hide rare aggregates.
- F. Take a look at info to make sure the number of counts is within the optical range for your sample (For Calibration standards, Regular image size: 1500-3000/ Large image size: 5000-9000).

- G. Create a calibration curve.
 - 1. Choose your mass calibration standard file(s).
 - 2. Clear the fitting curve which software automatically done on the histogram.
 - 3. Double click on each population to create a Gaussian fit (or click and drag over a peak).
 - 4. In Mass Calibration Section (bottom left), click "Create".
 - 5. Select the Calibration standard file(s) and enter the mass of your calibrants.
 - 6. Ensure the calibration curve has an R² of 1.000, and Max. mass error < 5% (recommend < 2%)
 - 7. Saving the calibration file and it will be added to Mass Calibrations Section.



- H. Apply the calibration curve to selected data files by double click the file of interest and then double click the calibration curve file in Mass Calibrations Section.
- I. Create Figures in DiscoverMP.
 - 1. Click the "Figures" tab and select "Add Figure"
 - 2. Choose 2D histogram, 3D histogram or Vertical histogram series.
 - 3. Drag the measurement file(s) and drop it in the figure set-up box.
 - 4. Deselect any aspects of the data you don't wish to include.
 - 5. Use the formatting option to change the graph title, axes, units and style.
- J. Save the workspace to preserve your analysis. [File -> Save Workspace]
- K. Export your figure(s).
- V. Shutdown
 - A. Close the AquireMP by clicking [File -> Quit]. The stage will return to the default position.
 - B. Turn off the electronics unit and vibration isolation bench (Press "Isolation" and then "Power").
 - C. Clean the objective with Lens Cleaning Tissue (Not Kimwipe!) and isopropanol when data collection is complete.
 - D. Close the lid to avoid dust accumulation.
 - E. Delete all your files in PC-1 and shutdown PC-1 after transfer them into PC-2
 - F. Shut down PC-2 when you finish your analysis.
 - G. Empty the water bucket inside the dehumidifier